


Intelligent

?LR@5MRGA



Classification & Composition

Classification method of intelligent instantaneous sterilizing equipment



Aj_qqj_a_rgmI_aampbjle rm
bj_cpcIrfc_rgle k crfmbq



Bjpcar_fc_rgle



Slbjpcar_fc_rgle



Aj_qqj_a_rgmI_aampbjle rm
bj_cpcIrbjgldcargmI k crfmbq

Qrc_k_bjgqldcargmI
kqrc_k_bjgqldcargmI'

Qsncpfc_rcb_u_rcp bjgqldca+
rgmI ku_rcp bjgqldcargmI'

Intelligent instantaneous sterilizing equipment is divided into 4 series

- +
- +
- +
- +



Rwnc 6gIrcjjjecIrgIqr_lrcIcmsq qrcpjpxgle cosjn kclr+ bjpcar
fc_rgle_l b qrc_k_bjgqldcargmI

Rwnc 66gIrcjjjecIrgIqr_lrcIcmsq qrcpjpxgle cosjn kclr+
bjpcar_fc_rgle_l b u_rcp bjgqldcargmI

Rwnc 666gIrcjjjecIrgIqr_lrcIcmsq qrcpjpxgle cosjn kclr+
glbjpcar_fc_rgle_l b qrc_k_bjgqldcargmI

Rwnc 6TgIrcjjjecIrgIqr_lrcIcmsq qrcpjpxgle cosjn kclr+
glbjpcar_fc_rgle_l b u_rcp bjgqldcargmI

K_gIrf_gIcp



Q5+/0...NJA

/...kq

MB@A-QQJ
MNA BN KMB@SQ
QCPTCP AJXCLR
BAQ

